



Results/Run 2/ Plate 1 (dmso + library at 0.2 μ M) QC Control wells




Zfactor 0.34994241




Zfactor 0.07936598

16

Results/Run 2/ Plate 2 (dmso + library at 2 μ M) QC Control wells




Zfactor 0.198704529




Zfactor -0.08355616

17

Results/Run 2/ Plate 3 (CBD + library at 0.2 μ M) QC Control wells




Zfactor 0.51979796




Zfactor -0.0396623

18

Results/Run 2/ Plate 4 (CBD + library at 2 μ M) QC Control wells



Zfactor 0.39309778




Zfactor 0.24674941


19

Results/Run 2/ Plate 5(Unperturbed) QC Control wells

1st plate to add CTB



AVG	430677.125
STD DEV	34643.7958
%CV	8.04402969




AVG	703781.427
STD DEV	93563.0908
%CV	13.2943393


20

Results/Run 2/ Plate 9 (Unperturbed) QC Control wells

Last plate to add CTB




AVG	437723.19
STD DEV	30927.8108
%CV	7.06560939




AVG	581847.201
STD DEV	95858.1036
%CV	16.4747899

21



Row/column variation on unperturbed plates from CTB
readout_plate 5

22



Row/column variation on unperturbed plates from CTB
readout_plate 6

23

Conclusions

- Z score_variation= we have high variability within the DMSO wells and staurosporine. Possible reasons for this variability:
 - Performance of wellmate = liquid being added on the wall/edges of the wells
 - New 384-well plate = we used Corning 3712 for run 1, and Corning 3764 for run 2 (3764 is supposed to be the replacement for 3712, they have overall the same characteristics, but they look quite different by “eye”).
- Z scores_overall= the reason why the separation is not large might be because Cal27 is not as sensitive to staurosporine, as Fadu is. Furthermore, we have drugs within the library that are consistently killing more cells than staurosporine (Bortezomib and dinaciclib).
- Wellmate= the possibility for most of the variability coming from the wellmate are also due:
 - Unperturbed plates showed a clear row-to-row variation
 - Because I noticed liquid in the wall/edges/top of the well after seeding the cells and adding the CTB, I had to spin the plates. This spinning could possibly explain the higher variability in the outer wells
 - To possibly address this, Margret will replace the tubing cassette for the wellmate

24

New 384-well plates

	Corning® 3712 384-Well – RUN 1	Corning® 3764 384-Well – RUN 2
Well Bottom	Flat and clear	Flat and clear
Well Volume	112µL	112µL
Recommended Working Volume	20 to 80µL	20 to 80µL
Cell Growth Area	0.06cm ²	0.06cm ²
Required pin tool adjustment?	N/A	No
Required wellmate adjustment?	N/A	YES

Link for the plates

<https://www.capitolscientific.com/Corning-3712-384-Well-x-112L-Clear-Flat-Bottom-Cell-Culture-Microplates-with-Lid-TC-Treated-BL#:~:text=This%20Corning%C2%AE%20384%2Dwell,volume%20of%20to%2050%C2%B5L>.

<https://ecatalog.corning.com/life-sciences/b2c/US/en/Microplates/Assay-Microplates/384-Well-Microplates/Corning%C2%AE-384-well-Black-Clear-and-White-Clear-Bottom-Polystyrene-Microplates/p/3764>

25